

## ARCHITECT INSTRUMENT SETTING

Apolipoprotein B - ApoB (BSBE)

Method : Immunoturbidimetric

<input checked="" type="radio"/> Genral	<input type="radio"/> Calibration	<input type="radio"/> Smartwash	<input type="radio"/> Result	<input type="radio"/> Interpretation	
Assay :	<input type="text" value="APOB"/>	Type :	<input type="text" value="Photometric"/>	Version :	<input type="text" value="1"/>
Number :	<input type="text"/>	Assay :	<input type="text" value="Enabled"/>	Date :	<input type="text"/>
				Time :	<input type="text"/>
				Operator :	<input type="text" value="ADMIN"/>

<input checked="" type="radio"/> Reaction Defination	<input type="radio"/> Reagent / Sample	<input type="radio"/> Vallidity Check				
Reaction mode :	<input type="text" value="End up"/>	Read Time				
Wavelength	<table style="width: 100%; border: none;"> <tr> <td style="text-align: center; border: none;">Primary</td> <td style="text-align: center; border: none;">Secondary</td> </tr> <tr> <td style="border: 1px solid black; text-align: center;">340/700</td> <td style="border: none;"></td> </tr> </table>	Primary	Secondary	340/700		Main : <input type="text" value="33 - 33"/>
Primary	Secondary					
340/700						
Last Required read	<input type="text" value="33"/>	Flex : <input type="text" value="0-0"/>				
Absorbance range	<input type="text" value="0.0000 - 3.0000"/>	Color correction : <input type="text" value="0-0"/>				
Sample blank type	<input type="text" value="Self"/>	Blank : <input type="text" value="10 - 10"/>				

<input type="radio"/> Reaction Defination	<input checked="" type="radio"/> Reagent / Sample	<input type="radio"/> Vallidity Check					
Reagent :	<input type="text" value="APOB"/>	Reagent Volume :	<table style="width: 100%; border: none;"> <tr> <td style="text-align: center; border: none;">R1</td> <td style="text-align: center; border: none;">R2</td> </tr> <tr> <td style="border: 1px solid black; text-align: center;">175</td> <td style="border: 1px solid black; text-align: center;">58</td> </tr> </table>	R1	R2	175	58
R1	R2						
175	58						
Diluent	<input type="text" value="Saline"/>	Water volume :	<input type="text" value="0"/>				
Diluent dispense mode :	<input type="text" value="Type0"/>	Dispense mode :	<table style="width: 100%; border: none;"> <tr> <td style="text-align: center; border: none;">Type0</td> <td style="text-align: center; border: none;">Type0</td> </tr> </table>	Type0	Type0		
Type0	Type0						
Dilution name	Sample	Diluted Sample	Diluent	Water	Dilution Factor	Default dilution	
STANDARD	<input type="text" value="1.6"/>	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="1:0"/>	<input checked="" type="radio"/>	
1:4	<input type="text" value="25"/>	<input type="text" value="2"/>	<input type="text" value="75"/>	<input type="text" value="0"/>	<input type="text" value="1:4"/>	<input type="radio"/>	

<input type="radio"/> Reaction Defination	<input type="radio"/> Reagent / Sample	<input checked="" type="radio"/> Vallidity Check	
Reaction check type :	<input type="text" value="None"/>	Read time B range :	<input type="text"/>
Read time A range :	<input type="text"/>	Minimum absorbance :	<input type="text"/>
Calculation limit :	<input type="text"/>	Rate linearity % :	<input type="text" value="0"/>

