

ARCHITECT INSTRUMENT SETTING

Glucose-6-Phosphate Dehydrogenase (Medicon)

Method : UV Kinetics

<input checked="" type="radio"/> Genral	<input type="radio"/> Calibration	<input type="radio"/> Smartwash	<input type="radio"/> Result	<input type="radio"/> Interpretation	
Assay :	<input type="text" value="G6PD"/>	Type :	<input type="text" value="Photometric"/>	Version :	<input type="text" value="1"/>
Number :	<input type="text"/>	Assay :	<input type="text" value="Enabled"/>	Date :	<input type="text"/>
				Time :	
				Operator :	<input type="text" value="ADMIN"/>

<input checked="" type="radio"/> Reaction Defination	<input type="radio"/> Reagent / Sample	<input type="radio"/> Vallidity Check
Reaction mode :	<input type="text" value="Rate Up"/>	
		Read Time
Wavelength	<input type="text" value="340/380"/> <div style="display: flex; justify-content: space-around; font-size: small;"> Primary Secondary </div>	Main : <input type="text" value="20-30"/>
Last Required read	<input type="text" value="30"/>	Flex : <input type="text" value="0-0"/>
Absorbance range	<input type="text" value="0.0000 - 3.0000"/>	Color correction : <input type="text" value="0-0"/>
Sample blank type	<input type="text" value="None"/>	Blank : <input type="text"/>

<input type="radio"/> Reaction Defination	<input checked="" type="radio"/> Reagent / Sample	<input type="radio"/> Vallidity Check					
Reagent :	<input type="text" value="G6PD"/>	Reagent Volume :	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 50%; text-align: center;">R1</td> <td style="width: 50%; text-align: center;">R2</td> </tr> <tr> <td style="text-align: center;">240</td> <td style="text-align: center;">20</td> </tr> </table>	R1	R2	240	20
R1	R2						
240	20						
Diluent	<input type="text" value="Saline"/>	Water volume :	<input type="text" value="0"/>				
Diluent dispense mode :	<input type="text" value="Type0"/>	Dispense mode :	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 50%; text-align: center;">Type0</td> <td style="width: 50%; text-align: center;">Type0</td> </tr> </table>	Type0	Type0		
Type0	Type0						
Dilution name	Sample	Diluted Sample	Diluent	Water	Dilution Factor	Default dilution	
STANDARD :	<input type="text" value="9.6"/>	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="1:1"/>	<input checked="" type="radio"/>	
	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="radio"/>	

<input type="radio"/> Reaction Defination	<input type="radio"/> Reagent / Sample	<input checked="" type="radio"/> Vallidity Check	
Reaction check type :	<input type="text" value="None"/>	Read time B range :	<input type="text"/>
Read time A range :	<input type="text"/>	Minimum absorbance :	<input type="text"/>
Calculation limit :	<input type="text"/>	Rate linearity % :	<input type="text" value="0"/>

<input type="radio"/> Genral	<input checked="" type="radio"/> Calibration	<input type="radio"/> Smartwash	<input type="radio"/> Result	<input type="radio"/> Interpretation	
Assay :	<input type="text" value="G6PD"/>	Assay Numbe :	<input type="text"/>	Date :	<input type="text"/>
Calibration method :	<input type="text" value="Factor"/>	Time :	<input type="text"/>	Operator :	<input type="text" value="ADMIN"/>
Factor :	<input type="text" value="4226.0"/>				

<input checked="" type="radio"/> Calibrators	<input type="radio"/> Volumes	<input type="radio"/> Intervals	<input type="radio"/> Validity checks		
Calibrator Set :	<input type="text" value="ADA"/>	Calibrator level :	<input type="text" value="Water"/>	Concentration :	<input type="text" value="0.00"/>
Blank :	<input type="text"/>	Cal 1 :	<input type="text"/>		
Replicates :	<input type="text" value="3"/>				

<input type="radio"/> Calibrators	<input checked="" type="radio"/> Volumes	<input type="radio"/> Intervals	<input type="radio"/> Validity checks		
Calibrator :	<input type="text" value="G6PD"/>				
Blank	<input type="text" value="Water"/>	<input type="text" value="3.5"/>	<input type="text" value="0.0"/>	<input type="text" value="0"/>	<input type="text" value="0"/>
Cal 1	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

<input type="radio"/> Genral	<input type="radio"/> Calibration	<input checked="" type="radio"/> Smartwash	<input type="radio"/> Result	<input type="radio"/> Interpretation	
Componentt	Reagent / Assay	Wash	Volume	Replicates	Wash Protocol
Sample Prob		Water			Madium Wash
Cuvette	TRIG	Detergent A	345		
Reagent Probe (R1)	All	Detergent A	345	1	
Reagent Probe (R1)	All	0.5% Acid Wash	345	1	
Reagent Probe (R1)	All	Water	345	2	
Reagent Probe (R2)	All	Detergent A	345	1	
Reagent Probe (R2)	All	0.5% Acid Wash	345	1	
Reagent Probe (R2)	All	Water	345	2	

<input type="radio"/> Genral	<input type="radio"/> Calibration	<input type="radio"/> Smartwash	<input checked="" type="radio"/> Result	<input type="radio"/> Interpretation
Linearty range :		0.00 - 3000.00		
Flag range specifaation :				

Result units				
Result Concentration unit :	<input type="text" value="U/L"/>	Correlation factor :	<input type="text" value="1.0000"/>	
Result decimal places :	<input type="text" value="2"/>	Intercept :	<input type="text" value="0.0000"/>	