

ARCHITECT INSTRUMENT SETTING

Lipoprotein (a) - Lp(a) (BSBE)
Method : Immunoturbidimetric

 Genral Calibration Smartwash Result InterpretationAssay : Type : Version : Number : Assay : Date : Time : Operator : Reaction Defination Reagent / Sample Vallidity CheckReaction mode :

Read Time

Wavelength Primary SecondaryMain : Last Required read Flex : Absorbance range Color correction : Sample blank type Blank : Reaction Defination Reagent / Sample Vallidity CheckReagent :

Reagent Volume :

R1

R2

Diluent

Water volume :

Diluent dispense mode :

Dispense mode :

Dilution name	Sample	Diluted Sample	Diluent	Water	Dilution Factor	Default dilution
STANDARD	1.5	0	0	0	1:0	<input checked="" type="radio"/>
1:10	10.0	2.0	90	0	1:10	<input type="radio"/>

 Reaction Defination Reagent / Sample Vallidity CheckReaction check type : Read time B range : Read time A range : Minimum absorbance : Calculation limit :

Rate linearity % :

<input type="radio"/> Genral	<input checked="" type="radio"/> Calibration	<input type="radio"/> Smartwash	<input type="radio"/> Result	<input type="radio"/> Interpretation
------------------------------	--	---------------------------------	------------------------------	--------------------------------------

Assay : Assay Numbe : Date :
 Calibration method : Operator :

<input checked="" type="radio"/> Calibrators	<input type="radio"/> Volumes	<input type="radio"/> Intervals	<input type="radio"/> Validity checks
--	-------------------------------	---------------------------------	---------------------------------------

Calibrator Set : Calibrator level : Concentration :
 Blank : Cal 1 :
 Replicates :

<input type="radio"/> Calibrators	<input checked="" type="radio"/> Volumes	<input type="radio"/> Intervals	<input type="radio"/> Validity checks
-----------------------------------	--	---------------------------------	---------------------------------------

Calibrator :

	Calibrator Level	Sample	Diluted sample	Diluent	Water
Blank	Water	1.5	0.0	0	0
Lp(a) CAL1	CAL 1	1.5	0.0	0	0
Lp(a) CAL2	CAL 2	1.5	0.0	0	0
Lp(a) CAL3	CAL 3	1.5	0.0	0	0
Lp(a) CAL4	CAL 4	1.5	0.0	0	0
Lp(a) CAL5	CAL 5	1.5	0.0	0	0

<input type="radio"/> Genral	<input type="radio"/> Calibration	<input checked="" type="radio"/> Smartwash	<input type="radio"/> Result	<input type="radio"/> Interpretation
------------------------------	-----------------------------------	--	------------------------------	--------------------------------------

<input type="radio"/> Genral	<input type="radio"/> Calibration	<input type="radio"/> Smartwash	<input checked="" type="radio"/> Result	<input type="radio"/> Interpretation
------------------------------	-----------------------------------	---------------------------------	---	--------------------------------------

Linearty range : 0.00 - 100.00
 Flag range specificaation :

Gender	Age	Normal Range	External range
Either	18 - 60 Years	0.00 - 30.00	

Result units

Result Concentration unit : Correlation factor :
 Result decimal places : Intercept :