

ARCHITECT INSTRUMENT SETTING

URINARY PROTEIN/CSF-UPRO (Medicon)

Method : Pyrogallol Red

<input checked="" type="radio"/> Genral	<input type="radio"/> Calibration	<input type="radio"/> Smartwash	<input type="radio"/> Result	<input type="radio"/> Interpretation	
Assay :	<input type="text" value="UPRO"/>	Type :	<input type="text" value="Photometric"/>	Version :	<input type="text" value="1"/>
Number :	<input type="text"/>	Assay :	<input type="text" value="Enabled"/>	Date :	<input type="text"/>
				Time :	<input type="text"/>
				Operator :	<input type="text" value="ADMIN"/>

<input checked="" type="radio"/> Reaction Defination	<input type="radio"/> Reagent / Sample	<input type="radio"/> Vallidity Check				
Reaction mode :	<input type="text" value="End up"/>	Read Time				
Wavelength	<table style="width: 100%; border: none;"> <tr> <td style="text-align: center;">Primary</td> <td style="text-align: center;">Secondary</td> </tr> <tr> <td style="border: 1px solid black; text-align: center;">604 / 804</td> <td style="border: 1px solid black;"></td> </tr> </table>	Primary	Secondary	604 / 804		Main : <input type="text" value="7 - 8"/>
Primary	Secondary					
604 / 804						
Last Required read	<input type="text" value="8"/>	Flex : <input type="text" value="0-0"/>				
Absorbance range	<input type="text" value="0.0000 - 3.0000"/>	Color correction : <input type="text" value="0-0"/>				
Sample blank type	<input type="text" value="None"/>	Blank : <input type="text"/>				

<input type="radio"/> Reaction Defination	<input checked="" type="radio"/> Reagent / Sample	<input type="radio"/> Vallidity Check					
Reagent :	<input type="text" value="UPRO"/>	Reagent Volume :	<table style="width: 100%; border: none;"> <tr> <td style="text-align: center;">R1</td> <td style="text-align: center;">R2</td> </tr> <tr> <td style="border: 1px solid black; text-align: center;">190</td> <td style="border: 1px solid black;"></td> </tr> </table>	R1	R2	190	
R1	R2						
190							
Diluent	<input type="text" value="Saline"/>	Water volume :	<input type="text" value="0"/>				
Diluent dispense mode :	<input type="text" value="Type0"/>	Dispense mode :	<input type="text" value="Type0"/>				
Dilution name	Sample	Diluted Sample	Diluent	Water	Dilution Factor	Default dilution	
STANDARD	<input type="text" value="2"/>	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="1:0"/>	<input checked="" type="radio"/>	
1:2	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="radio"/>	

<input type="radio"/> Reaction Defination	<input type="radio"/> Reagent / Sample	<input checked="" type="radio"/> Vallidity Check	
Reaction check type :	<input type="text" value="None"/>	Read time B range :	<input type="text"/>
Read time A range :	<input type="text"/>	Minimum absorbance :	<input type="text"/>
Calculation limit :	<input type="text"/>	Rate linearity % :	<input type="text" value="0"/>

<input type="radio"/> Genral	<input checked="" type="radio"/> Calibration	<input type="radio"/> Smartwash	<input type="radio"/> Result	<input type="radio"/> Interpretation	
Assay :	<input type="text" value="UPRO"/>	Assay Numbe :	<input type="text"/>	Date :	<input type="text"/>
Calibration method :	<input type="text" value="Linear"/>	Time :	<input type="text"/>	Operator :	<input type="text" value="ADMIN"/>

<input checked="" type="radio"/> Calibrators	<input type="radio"/> Volumes	<input type="radio"/> Intervals	<input type="radio"/> Validity checks
Calibrator Set :	<input type="text" value="UPRO"/>	Calibrator level :	Concentration :
		Blank :	<input type="text" value="0.00"/>
		Cal 1 :	<input type="text" value="UPRO CAL 1"/>
Replicates :	<input type="text" value="2"/>		

<input type="radio"/> Calibrators	<input checked="" type="radio"/> Volumes	<input type="radio"/> Intervals	<input type="radio"/> Validity checks		
Calibrator :	<input type="text" value="UPRO Cal"/>				
Blank	Water	17	0.0	0	0
Cal 1	UPRO Cal	17	0.0	0	0

<input type="radio"/> Genral	<input type="radio"/> Calibration	<input checked="" type="radio"/> Smartwash	<input type="radio"/> Result	<input type="radio"/> Interpretation	
Componentt	Reagent / Assay	Wash	Volume	Replicates	Wash Protocol
Sample Prob		Water			Madium Wash
Reagent Probe (R1)	All	Water	345	2	

<input type="radio"/> Genral	<input type="radio"/> Calibration	<input type="radio"/> Smartwash	<input checked="" type="radio"/> Result	<input type="radio"/> Interpretation
Linearty range :		0.75 - 200.00		
Flag range specifaaction :				
Gender	Age	Normal Range	External range	

Result units	
Result Concentration unit :	<input type="text" value="mg/L"/>
Result decimal places :	<input type="text" value="1"/>
Correlation factor :	<input type="text" value="1.0000"/>
Intercept :	<input type="text" value="0.0000"/>