

Roche Cobas 501/ Cobas 311 Instrument Setting

ACE (BSBE)

Method - Colorimetric

Analyze	Calib	Range	Other
Assay / Time / Point (C501)	Rate A	10	22 70 0 0
Assay / Time / Point (C311)	Rate A	10	15 50 0 0
Wavelength (2nd/Pri)	<input type="text" value="Cancle"/>	<input type="text" value="340"/>	
Sample Volume			Cassette Configuration
Norm.	<input type="text" value="14"/> <input type="text" value="0.0"/> <input type="text" value="0"/>		Code
Dec.	<input type="text" value="7"/> <input type="text" value="0.0"/> <input type="text" value="0"/>		Expiration Days <input type="text" value="99"/>
Inc.	<input type="text" value="28"/> <input type="text" value="0.0"/> <input type="text" value="0"/>		Reagent Volume
Dilution			R1 <input type="text" value="120"/> <input type="text" value="0"/> Inactive
<input checked="" type="radio"/> Water			R2 <input type="text" value="0"/> <input type="text" value="0"/>
<input type="radio"/> Diluent	<input type="text" value="0"/> <input type="text" value="1"/>		R3 <input type="text" value="0"/> <input type="text" value="0"/>
Linearity Limit	<input type="text" value="15"/> % <input type="text" value=""/>	% <input type="text" value=""/>	Bottle Setting
Prozone Limit	<input type="text" value="0"/> <input type="text" value="0"/>	<input type="text" value="0"/> <input type="text" value="0"/>	<input type="text" value="0"/> <input type="text" value="0"/> Inside <input type="text" value="0"/> <input type="text" value="0"/>
Abs. Limit	<input type="text" value="32000"/>	Decrease	
Cell Detergent	<input type="text" value="Detergent 1"/>	Stirring Level	<input type="text" value="1"/>
Stirring Setting		M1 M2 M3	SAVE
UP	<input type="text" value="Stirring"/> LOW	<input type="text" value="Stirring"/> <input type="text" value="Stirring"/> <input type="text" value="Stirring"/>	

Analyze	Calib	Range	Other
Calibration Type	<input type="text" value="Linear"/>		Auto Calibration
Point	<input type="text" value="2"/>		<input checked="" type="radio"/> Timeout
Span	<input type="text" value="2"/>		Cassette <input type="text" value="Cancle"/>
Weight	<input type="text" value="0"/>		<input type="text" value="0"/> Day
Update Type	<input type="text" value="None"/> <input type="text" value="0"/> <input type="text" value="0"/>		Changeover
			Cassette <input type="text" value="Cancle"/>
SD limit	<input type="text" value="999"/>		<input type="radio"/> QC Violation
Duplicate Limit	<input type="text" value="99"/> <input type="text" value="32000"/> Abs.		Method <input type="text" value="Blank"/>
Sensitivity Limit	<input type="text" value="-99999"/> <input type="text" value="99999"/>		Rule <input type="text" value="1 s"/>
S1 Abs. Limit	<input type="text" value="-32000"/> <input type="text" value="32000"/>		Control 1 <input type="text" value="None"/>
<input type="checkbox"/> Auto Masking			Control 2 <input type="text" value="None"/>
			Control 3 <input type="text" value="None"/>
			SAVE

Analyze	Calib	Range	Other
---------	-------	-------	-------

Application Code

Unit U/L

Report Name ACE

Data Mode Active

Automatic Rerun

Technical Kimit

5	150
---	-----

Repeat Limit

-99999	99999
--------	-------

Control Interval Time

Automatic QC On Board Stability

Qualitative

(1)	0		
(2)	0		
(3)	0		
(4)	0		
(5)	0		
(6)	0		

L 0

H 0

I 0

Default

Sex
 Male Female

Range
 Range 1 Range 2 Range 3

SAVE

Analyze	Calib	Range	Other
---------	-------	-------	-------

	(1)	(2)	(3)	(4)	(5)	(6)
Calibrator Code	<input style="width: 50px;" type="text"/>	<input style="width: 50px;" type="text"/>	<input style="width: 50px;" type="text"/>	<input style="width: 50px;" type="text"/>	<input style="width: 50px;" type="text"/>	<input style="width: 50px;" type="text"/>
Concentration	0.00	Cal 1				

Rack No. - Pos.

Sample Volume	14.0	14.0				
Diluted S. Volume	0.0	0.0	0.0	0.0	0.0	0.0
Diluent Volume	0	0	0	0	0	0

SAVE

Bottle Setting

Cassette Type A

Bottle

a	R2		
b	R1		
c	R1		

Cancel

OK