

Roche Cobas 501/ Cobas 311 Instrument Setting

Serum Copper (BSBE)

Method - Colorimetric

Analyze	Calib	Range	Other									
Assay / Time / Point (C501)	2 Point End	10	30									
Assay / Time / Point (C311)	2 Point End	10	20									
Wavelength (2nd/Pri)	700	600										
Sample Volume	Cassette Configuration											
Norm. <table border="1" style="display: inline-table; border-collapse: collapse;"><tr><td>7</td><td>0.0</td><td>0</td></tr><tr><td>3.5</td><td>0.0</td><td>0</td></tr><tr><td>14</td><td>0.0</td><td>0</td></tr></table>	7	0.0	0	3.5	0.0	0	14	0.0	0	Code		
7	0.0	0										
3.5	0.0	0										
14	0.0	0										
Dec.	Expiration Days <input style="width: 50px;" type="text" value="99"/>											
Inc.	Reagent Volume											
Dilution	R1 <table border="1" style="display: inline-table; border-collapse: collapse;"><tr><td>105</td><td>0</td><td>Inactive</td></tr><tr><td>0</td><td>0</td><td></td></tr><tr><td>35</td><td>0</td><td></td></tr></table>		105	0	Inactive	0	0		35	0		Bottle Setting
105	0	Inactive										
0	0											
35	0											
<input checked="" type="radio"/> Water	R2											
<input type="radio"/> Diluent	R3											
Linearity Limit	<input style="width: 50px;" type="text" value="15"/> %	<input style="width: 50px;" type="text" value="0"/> %	<input style="width: 50px;" type="text" value="0"/> %									
Prozone Limit	<input style="width: 50px;" type="text" value="0"/> <input style="width: 50px;" type="text" value="0"/>	<input style="width: 50px;" type="text" value="0"/> <input style="width: 50px;" type="text" value="0"/>	<input style="width: 50px;" type="text" value="0"/> <input style="width: 50px;" type="text" value="0"/> Inside <input style="width: 50px;" type="text" value="0"/> <input style="width: 50px;" type="text" value="0"/>									
Abs. Limit	<input style="width: 50px;" type="text" value="32000"/>	Increase										
Cell Detergent	<input style="width: 100px;" type="text" value="Detergent 1"/>		Stirring Level <input style="width: 50px;" type="text" value="1"/>									
Stirring Setting	M1 M2 M3											
UP <input style="width: 50px;" type="text" value="Stirring"/>	LOW <input style="width: 50px;" type="text" value="Stirring"/>	<input style="width: 50px;" type="text" value="Stirring"/>	<input style="width: 50px;" type="text" value="Stirring"/>									
			SAVE									

Analyze	Calib	Range	Other
Calibration Type	Linear		
Point	<input style="width: 50px;" type="text" value="2"/>		
Span	<input style="width: 50px;" type="text" value="2"/>		
Weight	<input style="width: 50px;" type="text" value="0"/>		
Update Type	None	<input style="width: 50px;" type="text" value="0"/>	<input style="width: 50px;" type="text" value="0"/>
SD limit	<input style="width: 50px;" type="text" value="999"/>		
Duplicate Limit	<input style="width: 50px;" type="text" value="99"/>	<input style="width: 50px;" type="text" value="32000"/>	Abs.
Sensitivity Limit	<input style="width: 50px;" type="text" value="-99999"/>	<input style="width: 50px;" type="text" value="99999"/>	
S1 Abs. Limit	<input style="width: 50px;" type="text" value="-32000"/>	<input style="width: 50px;" type="text" value="32000"/>	
<input type="checkbox"/> Auto Masking			
Auto Calibration	<input checked="" type="radio"/> Timeout		
Cassette	<input style="width: 50px;" type="text" value="Cancle"/>		
	<input style="width: 50px;" type="text" value="0"/>	Day	
Changeover	<input style="width: 50px;" type="text" value="Cancle"/>		
Cassette	<input style="width: 50px;" type="text" value="Cancle"/>		
QC Violation	<input type="radio"/>		
Method	Blank		
Rule	1 s		
Control 1	None		
Control 2	None		
Control 3	None		
			SAVE

Analyze	Calib	Range	Other
---------	-------	-------	-------

Application Code

Unit

Report Name

Data Mode

Automatic Rerun

Technical Kimit	0.3	60
Repeat Limit	-99999	99999

Control Interval Time

Automatic QC On Board Stability

<input type="checkbox"/> Qualitative	
(1)	0
(2)	0
(3)	0
(4)	0
(5)	0
(6)	

L	0
H	0
I	0

Default

Sex
 Male Female

Range
 Range 1 Range 2 Range 3

Analyze	Calib	Range	Other
---------	-------	-------	-------

	(1)	(2)	(3)	(4)	(5)	(6)
Calibrator Code	<input style="width: 50px;" type="text"/>	<input style="width: 50px;" type="text"/>	<input style="width: 50px;" type="text"/>	<input style="width: 50px;" type="text"/>	<input style="width: 50px;" type="text"/>	<input style="width: 50px;" type="text"/>
Concentration	0.00	Cal 1				

Rack No. - Pos.

Sample Volume	7.0	7.0				
Diluted S. Volume	0.0	0.0	0.0	0.0	0.0	0.0
Diluent Volume	0	0	0	0	0	0

Bottle Setting

Cassette Type A

Bottle			
a	R2		
b	R1		
c	R1		

RedCells Diagnostics Pvt Ltd

www.redcellsdiagnostics.com

Feb-19