

Roche Cobas 501/ Cobas 311 Instrument Setting

Serum Zinc (BSBE)

Method - Colorimetric

Analyze	Calib	Range	Other
Assay / Time / Point (C501)	1 Point	10	48 0 0 0
Assay / Time / Point (C311)	1 Point	10	35 0 0 0
Wavelength (2nd/Pri)	700	570	
Sample Volume	Cassette Configuration		
Norm. 6 0.0 0	Code		
Dec. 3 0.0 0	Expiration Days 99		
Inc. 12 0.0 0	Reagent Volume		
Dilution			Bottle Setting
<input checked="" type="radio"/> Water			
<input type="radio"/> Diluent	0	1	
Linearity Limit	15 %	%	
Prozone Limit	0 0	0 0	0 0 Inside 0 0
Abs. Limit	32000	Increase	
Cell Detergent	Detergent 1		Stirring Level 1
Stirring Setting	M1 M2 M3		SAVE
UP	Stirring LOW	Stirring Stirring	

Analyze	Calib	Range	Other
Calibration Type	Linear		
Point	2		
Span	2		
Weight	0		
Update Type	None 0 0		
Auto Calibration	<input checked="" type="radio"/> Timeout		
Cassette	Cancle		0 Day
Changeover	Cassette Cancle		
SD limit	999		
Duplicate Limit	99	32000	Abs.
Sensitivity Limit	-99999	99999	
S1 Abs. Limit	-32000	32000	
<input type="checkbox"/> Auto Masking	<input type="radio"/> QC Violation		
	Method	Blank	
	Rule	1 s	
	Control 1	None	
	Control 2	None	
	Control 3	None	
	SAVE		

Analyze	Calib	Range	Other																								
Application Code <input style="width: 60px;" type="text"/> Unit $\mu\text{g/dL}$ Report Name <input style="width: 100px;" type="text" value="Serum Zinc"/> Data Mode <input style="width: 60px;" type="text" value="Active"/>		Expected Values Male <table border="1" style="width: 100%; border-collapse: collapse;"> <tr><td style="width: 25%;">99</td><td style="width: 25%;">Year</td><td style="width: 25%;">70</td><td style="width: 25%;">115</td></tr> <tr><td>100</td><td>Year</td><td>-99999</td><td>999999</td></tr> <tr><td></td><td></td><td>-99999</td><td>999999</td></tr> </table> Female <table border="1" style="width: 100%; border-collapse: collapse;"> <tr><td style="width: 25%;">99</td><td style="width: 25%;">Year</td><td style="width: 25%;">70</td><td style="width: 25%;">115</td></tr> <tr><td>100</td><td>Year</td><td>-99999</td><td>999999</td></tr> <tr><td></td><td></td><td>-99999</td><td>999999</td></tr> </table>		99	Year	70	115	100	Year	-99999	999999			-99999	999999	99	Year	70	115	100	Year	-99999	999999			-99999	999999
99	Year	70	115																								
100	Year	-99999	999999																								
		-99999	999999																								
99	Year	70	115																								
100	Year	-99999	999999																								
		-99999	999999																								
<input type="checkbox"/> Automatic Rerun Technical Kimit <table border="1" style="width: 100%; border-collapse: collapse;"><tr><td style="width: 50%;">6.92</td><td style="width: 50%;">500</td></tr></table> Repeat Limit <table border="1" style="width: 100%; border-collapse: collapse;"><tr><td style="width: 50%;">-99999</td><td style="width: 50%;">99999</td></tr></table>		6.92	500	-99999	99999																						
6.92	500																										
-99999	99999																										
<input type="checkbox"/> Control Interval Time <input type="checkbox"/> Automatic QC On Board Stability																											
<input type="checkbox"/> Qualitative <table border="1" style="width: 100%; border-collapse: collapse;"> <tr><td style="width: 10%;">(1)</td><td style="width: 10%;">0</td><td style="width: 10%;">L</td><td style="width: 10%;">0</td></tr> <tr><td>(2)</td><td>0</td><td>H</td><td>0</td></tr> <tr><td>(3)</td><td>0</td><td>I</td><td>0</td></tr> <tr><td>(4)</td><td>0</td><td></td><td></td></tr> <tr><td>(5)</td><td>0</td><td></td><td></td></tr> <tr><td>(6)</td><td></td><td></td><td></td></tr> </table>		(1)	0	L	0	(2)	0	H	0	(3)	0	I	0	(4)	0			(5)	0			(6)				Default Sex <input checked="" type="radio"/> Male <input type="radio"/> Female Range <input checked="" type="radio"/> Range 1 <input type="radio"/> Range 2 <input type="radio"/> Range 3	
(1)	0	L	0																								
(2)	0	H	0																								
(3)	0	I	0																								
(4)	0																										
(5)	0																										
(6)																											
			SAVE																								

Analyze	Calib	Range	Other																					
<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 15%;"></td> <td style="width: 15%; text-align: center;">(1)</td> <td style="width: 15%; text-align: center;">(2)</td> <td style="width: 15%; text-align: center;">(3)</td> <td style="width: 15%; text-align: center;">(4)</td> <td style="width: 15%; text-align: center;">(5)</td> <td style="width: 15%; text-align: center;">(6)</td> </tr> <tr> <td>Calibrator Code</td> <td><input style="width: 60px;" type="text"/></td> <td><input style="width: 60px;" type="text"/></td> <td><input style="width: 60px;" type="text"/></td> <td><input style="width: 60px;" type="text"/></td> <td><input style="width: 60px;" type="text"/></td> <td><input style="width: 60px;" type="text"/></td> </tr> <tr> <td>Concentration</td> <td>0.00</td> <td>Cal 1</td> <td></td> <td></td> <td></td> <td></td> </tr> </table>					(1)	(2)	(3)	(4)	(5)	(6)	Calibrator Code	<input style="width: 60px;" type="text"/>	<input style="width: 60px;" type="text"/>	<input style="width: 60px;" type="text"/>	<input style="width: 60px;" type="text"/>	<input style="width: 60px;" type="text"/>	<input style="width: 60px;" type="text"/>	Concentration	0.00	Cal 1				
	(1)	(2)	(3)	(4)	(5)	(6)																		
Calibrator Code	<input style="width: 60px;" type="text"/>	<input style="width: 60px;" type="text"/>	<input style="width: 60px;" type="text"/>	<input style="width: 60px;" type="text"/>	<input style="width: 60px;" type="text"/>	<input style="width: 60px;" type="text"/>																		
Concentration	0.00	Cal 1																						
Rack No. - Pos. <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 15%;">Sample Volume</td> <td style="width: 15%;">7.0</td> <td style="width: 15%;">7.0</td> <td style="width: 15%;"></td> <td style="width: 15%;"></td> <td style="width: 15%;"></td> <td style="width: 15%;"></td> </tr> <tr> <td>Diluted S. Volume</td> <td>0.0</td> <td>0.0</td> <td>0.0</td> <td>0.0</td> <td>0.0</td> <td>0.0</td> </tr> <tr> <td>Diluent Volume</td> <td>0</td> <td>0</td> <td>0</td> <td>0</td> <td>0</td> <td>0</td> </tr> </table>				Sample Volume	7.0	7.0					Diluted S. Volume	0.0	0.0	0.0	0.0	0.0	0.0	Diluent Volume	0	0	0	0	0	0
Sample Volume	7.0	7.0																						
Diluted S. Volume	0.0	0.0	0.0	0.0	0.0	0.0																		
Diluent Volume	0	0	0	0	0	0																		
			SAVE																					

Bottle Setting																
Cassette Type A <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td colspan="4">Bottle</td> </tr> <tr> <td style="width: 10%;">a</td> <td style="width: 15%;">R2</td> <td style="width: 15%;"></td> <td style="width: 15%;"></td> </tr> <tr> <td>b</td> <td>R1</td> <td></td> <td></td> </tr> <tr> <td>c</td> <td>R1</td> <td></td> <td></td> </tr> </table>	Bottle				a	R2			b	R1			c	R1		
Bottle																
a	R2															
b	R1															
c	R1															
Cancel	OK															