

DADE DIMENSION PARAMETER SETTINGS

G-cell

ADA

Channel: *	Name:	XADA	Mode:	ABSORBENCE	Std. curve:	LINEAR
Delivery	Time	Component 1	Component 2	Component 3	Chase	Mix
R1	-60.0 sec	(A)150 µl	0 µl	0 µl	00 µl	NONE
S1	0.0 sec	4 µl	0 µl	0 µl	10 µl	MODERATE
R2	120 sec	(B)75 µl	0 µl	0 µl	20 µl	MODERATE
R3	0 sec	0 µl	0 µl	0 µl	NONE µl	
Photometry	Time	Cartridge configuration				
P1	300 sec	Component	(A)	(A)	(A)	(B) (B)
P2	420 sec	Number of tests	20	20	20	20 40 40
P3	sec	Well life hours	720	720	720	720 720 720
P4	sec	On board life 720 hours	Calibration 2160 hours			
Calculation	Mode	Measuring filter	546 nm	Blanking filter	800 nm	
P1	Time 300 sec	Dilution	1.0	IOD	0	(Inactive)
P2	Time 420 sec	Dilution	1.0	FOD	0	(Inactive)
	1	2	3	4	5	6
	*	*	*	*	*	*
Decimal	:	2				
Units	:	U/L				
Range	:	Serum 4-20 CSF 0-5 Pleural fluid 4-40				
ASSAY Range	:	0-176				