

DADE DIMENSION PARAMETER SETTINGS

G-cell

ALB

Channel:	*	Name:	XALB	Mode:	ABSORBENCE	Std. curve:	LINEAR
Delivery	Time	Component 1	Component 2	Component 3	Chase	Mix	
R1	-60.0 sec	(A)200 µl	0 µl	0 µl	20 µl	NONE	
S1	0.0 sec	2 µl	0 µl	0 µl	10 µl	MODERATE	
R2	sec	(B) µl	0 µl	0 µl	0 µl		
R3	0 sec	0 µl	0 µl	0 µl	NONE µl		
Photometry	Time	Cartridge configuration					
P1	-30 sec	Component	(A)	(A)	(A)	(A)	(A)
P2	80 sec	Number of tests	18	18	18	18	18
P3	sec	Well life hours	720	720	720	720	720
P4	sec	On board life	720 hours	Calibration	2160 hours		
Calculation							
Mode	END POINT	Measuring filter	600 nm	Blanking filter	700 nm		
P1	Time -30 sec	Dilution	-1.0	IOD	0	(Inactive)	
P2	Time 180 sec	Dilution	1.0	FOD	0	(Inactive)	
	1	2	3	4	5	6	
	*	*	*	*	*	*	
Decimal	2						
Units	g/L						
Interval	Male						
Range	3.8-4.4						