

DADE DIMENSION PARAMETER SETTINGS

G-cell

ALP

Channel: *	Name:	XALP	Mode:	ABSORBENCE	Std. curve:	LINEAR
Delivery	Time	Component 1	Component 2	Component 3	Chase	Mix
R1	-60.0 sec	(A) 160 µl	0 µl	0 µl	µl	NONE
S1	0.0 sec	3.0 µl	0 µl	0 µl	10 µl	MODERATE
R2	120 sec	(B) 40 µl	0 µl	0 µl	20 µl	MODERATE
R3	0 sec	0 µl	0 µl	0 µl	NONE µl	
Photometry	Time	Cartridge configuration				
P1	180 sec	Component	(A)	(A)	(A)	(B) (B)
P2	360 sec	Number of tests	15	15	15	30 30
P3	sec	Well life hours	720	720	720	720 720
P4	sec	On board life 720 hours	Calibration 2160 hours			
Calculation						
Mode	RATE	Measuring filter	405 nm	Blanking filter	600 nm	
P1	Time 180 sec	Dilution	1.0	IOD	0	(Inactive)
P2	Time 360 sec	Dilution	1.0	FOD	0	(Inactive)
1	2	3	4	5	6	
*	*	*	*	*	*	*
Decimal	: 2					
Units	: U/L					
Range	: 45 – 135					
ASSAY Range	: 0-1705					