

# DADE DIMENSION PARAMETER SETTINGS

# G-cell

## ALP

Channel: *	Name:	XALP	Mode:	ABSORBENCE	Std. curve:	LINEAR
<b>Delivery</b>	<b>Time</b>	<b>Component 1</b>	<b>Component 2</b>	<b>Component 3</b>	<b>Chase</b>	<b>Mix</b>
R1	-60.0 sec	(A) 160 µl	0 µl	0 µl	µl	NONE
S1	0.0 sec	3.0 µl	0 µl	0 µl	10 µl	MODERATE
R2	120 sec	(B) 40 µl	0 µl	0 µl	20 µl	MODERATE
R3	0 sec	0 µl	0 µl	0 µl	NONE µl	
<b>Photometry</b>	<b>Time</b>	<b>Cartridge configuration</b>				
P1	180 sec	Component	(A)	(A)	(A)	(B) (B)
P2	360 sec	Number of tests	15	15	15	15 30 30
P3	sec	Well life hours	720	720	720	720 720 720
P4	sec	On board life 720 hours	Calibration 2160 hours			
<b>Calculation</b>						
<b>Mode</b>	RATE	<b>Measuring filter</b>	405 nm	<b>Blanking filter</b>	600 nm	
P1	Time 180 sec	Dilution	1.0	IOD	0	(Inactive)
P2	Time 360 sec	Dilution	1.0	FOD	0	(Inactive)
	1	2	3	4	5	6
	*	*	*	*	*	*
<b>Decimal</b>	: 2					
<b>Units</b>	: U/L					
<b>Range</b>	: 45 – 135					
<b>ASSAY Range</b>	: 0-1705					