

DADE DIMENSION PARAMETER SETTINGS

G-cell

AMY

Channel: *	Name:	XAMY	Mode:	ABSORBENCE	Std. curve:	LINEAR
Delivery	Time	Component 1	Component 2	Component 3	Chase	Mix
R1	-60.0 sec	(A) 150 µl	0 µl	0 µl	20 µl	NONE
S1	0.0 sec	4 µl	0 µl	0 µl	10 µl	MODERATE
R2	120 sec	(B) 54 µl	0 µl	0 µl	20 µl	MODERATE
R3	0 sec	0 µl	0 µl	0 µl	NONE µl	
Photometry	Time	Cartridge configuration				
P1	180 sec	Component	(A)	(A)	(A)	(B)
P2	300 sec	Number of tests	12	12	12	48
P3	sec	Well life hours	720	720	720	720
P4	sec	On board life 720 hours	Calibration 2160 hours			
Calculation						
Mode	RATE	Measuring filter	405 nm	Blanking filter	700 nm	
P1	Time 180 sec	Dilution	1.0	IOD	0	(Inactive)
P2	Time 300 sec	Dilution	1.0	FOD	0	(Inactive)
	1	2	3	4	5	6
	*	*	*	*	*	*
Decimal	2					
Units	U/L					
Interval	25 – 104					
Range	3500					