

DADE DIMENSION PARAMETER SETTINGS

G-cell

APOA1

Channel:	*	Name:	XAPA	Mode:	TURBIDIMETR IC	Std. curve:	LOGIT
Delivery	Time	Component 1	Component 2	Component 3	Chase	Mix	
R1	-60.0 sec	(A)165 µl	0 µl	0 µl	20 µl	NONE	
S1	0.0 sec	1.5 µl	0 µl	0 µl	10 µl	MODERATE	
R2	100 sec	(B)55 µl	0 µl	0 µl	20 µl	MODERATE	
R3	0 sec	0 µl	0 µl	0 µl	0 µl		
Photometry	Time	Cartridge configuration					
P1	80 sec	Component	(A)	(A)	(A)	(A)	(B)
P2	400 sec	Number of tests	20	20	20	20	40
P3	sec	Well life hours	720	720	720	720	720
P4	sec	On board life 720 hours	Calibration 2160 hours				
Calculation							
Mode	END POINT	Measuring filter	340 nm	Blanking filter	700 nm		
P1	Time 80 sec	Dilution	-1.0	IOD	0	(Inactive)	
P2	Time 400 sec	Dilution	1.0	FOD	0	(Inactive)	
	1	2	3	4	5	6	
	*	*	*	*	*	*	
Decimal	2						
Units	mg/L						
Interval	120 – 176						
Assay Range	11-240						