

DADE DIMENSION PARAMETER SETTINGS

G-cell

CKMB

Channel: *	Name:	XCKB	Mode:	ABSORBENCE	Std. curve:	LINEAR
Delivery						
	Time	Component 1	Component 2	Component 3	Chase	Mix
R1	-60.0 sec	(A)160 µl	0 µl	0 µl	20 µl	NONE
S1	0.0 sec	8 µl	0 µl	0 µl	10 µl	MODERATE
R2	180 sec	(B)40 µl	0 µl	0 µl	20 µl	MODERATE
R3	0 sec	0 µl	0 µl	0 µl	NONE µl	
Photometry						
	Time	Cartridge configuration				
P1	300 sec	Component	(A)	(A)	(A)	(B) ()
P2	400 sec	Number of tests	15	15	15	60 0
P3	sec	Well life hours	720	720	720	720 720
P4	sec	On board life 720 hours	Calibration 2160 hours			
Calculation						
Mode	RATE	Measuring filter	340 nm	Blanking filter	405 nm	
P1	Time 300 sec	Dilution	1.0	IOD	0	(Inactive)
P2	Time 400 sec	Dilution	1.0	FOD	0	(Inactive)
	1	2	3	4	5	6
	*	*	*	*	*	*
Decimal	0					
Units	U/L					
Interval	0 – 25					
Assay Range	600					