

DADE DIMENSION PARAMETER SETTINGS

G-cell

CRE

Channel:	*	Name:	XCRE	Mode:	ABSORBENCE	Std. curve:	LINEAR
Delivery							
		Time	Component 1	Component 2	Component 3	Chase	Mix
R1		-60.0 sec	(A)160 µl	0 µl	0 µl	20 µl	NONE
S1		0.0 sec	8 µl	0 µl	0 µl	10 µl	MODERATE
R2		180 sec	(B)40 µl	0 µl	0 µl	20 µl	MODERATE
R3		0 sec	0 µl	0 µl	0 µl	NONE µl	
Photometry							
		Time	Cartridge configuration				
P1		240 sec	Component	(A)	(A)	(A)	(A) (B) (B)
P2		300 sec	Number of tests	15	15	15	15 30 30
P3		sec	Well life hours	720	720	720	720 720 720
P4		sec	On board life 720 hours	Calibration 2160 hours			
Calculation							
Mode		RATE	Measuring filter	510 nm	Blanking filter	660 nm	
P1		Time 240 sec	Dilution	1.0	IOD	0	(Inactive)
P2		Time 300 sec	Dilution	1.0	FOD	0	(Inactive)
		1	2	3	4	5	6
		*	*	*	*	*	*
Decimal		0					
Units		mg/dL					
Interval		0.8-1.4					
Assay Range		35					