

DADE DIMENSION PARAMETER SETTINGS



CRP

Channel:	*	Name:	XCRP	Mode:	ABSORBENC E	Std. curve:	Logit
Delivery							
		Time	Component 1		Component 2		Component 3
			Chase	Mix			
R1		-60.0 sec	(A)180 µl	0 µl	0 µl	20 µl	NONE
S1		0.0 sec	13 µl	0 µl	0 µl	10 µl	MODERATE
						20 µl	
R2		120 sec	(B)60 µl	0 µl	0 µl		MODERATE
R3		0 sec	0 µl	0 µl	0 µl		
Photometry							
		Time	Cartridge configuration				
P1		60 sec	Component	(A)	(A)	(A)	(A) (B) (B)
P2		420 sec	Number of tests	15	15	15	15 30 30
P3		sec	Well life hours	720	720	720	720 720 720
P4		sec	On board life 720 hours	Calibration 2160 hours			
Calculation							
Mode	END		Measuring filter	340 nm		Blanking filter	700 nm
P1	Time 60 sec		Dilution	-1.0		IOD	0 (Inactive)
P2	Time 420 sec		Dilution	1.0		FOD	0 (Inactive)
	1	2	3	4	5	6	
	*	*	*	*	*		
Decimal	2						
Units	mg/dL						
Interval	>0.8						
Range	205						