

DADE DIMENSION PARAMETER SETTINGS

G-cell

CYS-C

Cat. No: GB350S/GS351S

Channel: *	Name:	XCYS-C	Mode:	ABSORBENCE	Std. curve:	LINEAR
Delivery						
	Time	Component 1	Component 2	Component 3	Chase	Mix
R1	-60.0 sec	(A)150 µl	0 µl	0 µl	20 µl	MODERATE
S1	0.0 sec	2.5 µl	0 µl	0 µl	10 µl	MODERATE
R2	180 sec	(B)30 µl	0 µl	0 µl	20 µl	MODERATE
R3	0 sec	0 µl	0 µl	0 µl	NONE µl	
Photometry						
	Time	Cartridge configuration				
P1	220 sec	Component	(A)	(A)	(A)	(B) (B)
P2	450 sec	Number of tests	15	15	15	30 30
P3	sec	Well life hours	720	720	720	720 720
P4	sec	On board life 720 hours	Calibration 2160 hours			
Calculation						
Mode	END	Measuring filter	577 nm	Blanking filter	NONE nm	
P1	Time 240 sec	Dilution	-1.0	IOD	0	(Inactive)
P2	Time 450 sec	Dilution	1.0	FOD	0	(Inactive)
	1	2	3	4	5	6
	*	*	*			
Decimal	2					
Units	mg/L					
Interval	>1.5					
Range	0-8.00					