

DADE DIMENSION PARAMETER SETTINGS

G-cell

HCY

Channel:	*	Name:	XHCY	Mode:	ABSORBENCE	Std. curve:	LINEAR
Delivery	Time	Component 1	Component 2	Component 3	Chase	Mix	
R1	-60.0 sec	(A)150 µl	0 µl	0 µl	20 µl	NONE	
S1	0.0 sec	8 µl	0 µl	0 µl	10 µl	MODERATE	
R2	240 sec	(B)39 µl	0 µl	0 µl	20 µl	MODERATE	
R3	0 sec	0 µl	0 µl	0 µl	NONE µl		
Photometry	Time	Cartridge configuration					
P1	400 sec	Component	(A)	(A)	(A)	(A)	(B) (B)
P2	550 sec	Number of tests	15	15	15	15	30 30
P3	sec	Well life hours	720	720	720	720	720 720
P4	sec	On board life	720 hours	Calibration	2160 hours		
Calculation							
Mode	RATE	Measuring filter	340	nm	Blanking filter	700	nm
P1	Time 400 sec	Dilution	1.0		IOD	0	(Inactive)
P2	Time 550 sec	Dilution	1.0		FOD	0	(Inactive)
	1	2	3	4	5	6	
	*	*	*	*	*	*	*
Decimal	2						
Units	umol/L						
Interval	0-15						
Range	0-50						