

# DADE DIMENSION PARAMETER SETTINGS

# G-cell

## TBIL

Channel:	*	Name:	XTBIL	Mode:	ABSORBENCE	Std. curve:	LINEAR
<b>Delivery</b>							
	<b>Time</b>		<b>Component 1</b>	<b>Component 2</b>	<b>Component 3</b>	<b>Chase</b>	<b>Mix</b>
R1	-60.0 sec	(A)160 µl	0 µl	0 µl	0 µl	20 µl	NONE
S1	0.0 sec	8 µl	0 µl	0 µl	0 µl	10 µl	MODERATE
R2	120 sec	(B)40 µl	0 µl	0 µl	0 µl	20 µl	MODERATE
R3	0 sec	0 µl	0 µl	0 µl	0 µl	NONE µl	
<b>Photometry</b>							
	<b>Time</b>		<b>Cartridge configuration</b>				
P1	80 sec	Component	(A)	(A)	(A)	(A)	(B) (B)
P2	400 sec	Number of tests	15	15	15	15	30 30
P3	sec	Well life hours	720	720	720	720	720 720
P4	sec	On board life	720 hours	Calibration	2160 hours		
<b>Calculation</b>							
<b>Mode</b>	END POINT	<b>Measuring filter</b>	540 nm	<b>Blanking filter</b>	700 nm		
P1	Time 80 sec	Dilution	-1.0	IOD	0	(Inactive)	
P2	Time 400 sec	Dilution	1.0	FOD	0	(Inactive)	
	1	2	3	4	5	6	
	*	*	*	*	*	*	
<b>Decimal</b>	2						
<b>Units</b>	Mg/dL						
<b>Interval</b>	0.6-1.0						
<b>Assay Range</b>	30						