

# DADE DIMENSION PARAMETER SETTINGS

# G-cell

## TG

Channel: *	Name:	XTG	Mode:	ABSORBENCE	Std. curve:	LINEAR
<b>Delivery</b>	<b>Time</b>	<b>Component 1</b>	<b>Component 2</b>	<b>Component 3</b>	<b>Chase</b>	<b>Mix</b>
R1	-60.0 sec	(A)200 µl	0 µl	0 µl	20 µl	NONE
S1	0.0 sec	2 µl	0 µl	0 µl	20 µl	MODERATE
R2	sec	(B) µl	0 µl	0 µl	0 µl	
R3	0 sec	0 µl	0 µl	0 µl	NONE µl	
<b>Photometry</b>	<b>Time</b>	<b>Cartridge configuration</b>				
P1	-30 sec	Component	(A)	(A)	(A)	(A)
P2	300 sec	Number of tests	15	15	15	15
P3	sec	Well life hours	720	720	720	720
P4	sec	On board life 720 hours	Calibration 2160 hours			
<b>Calculation</b>	<b>Mode</b>	<b>END POINT</b>	<b>Measuring filter</b>	540 nm	<b>Blanking filter</b>	700 nm
P1	Time -30 sec	Dilution	-1.0	IOD	0	(Inactive)
P2	Time 300 sec	Dilution	1.0	FOD	0	(Inactive)
	1	2	3	4	5	6
	*	*	*	*	*	*
<b>Decimal</b>	2					
<b>Units</b>	Mg/dL					
<b>Interval</b>	60-150					
<b>Range</b>	1000					