

DADE DIMENSION PARAMETER SETTINGS

G-cell

TP

Channel: *	Name:	XTP	Mode:	ABSORBENCE	Std. curve:	LINEAR
Delivery	Time	Component 1	Component 2	Component 3	Chase	Mix
R1	-60.0 sec	(A)200 µl	0 µl	0 µl	20 µl	NONE
S1	0.0 sec	4 µl	0 µl	0 µl	20 µl	MODERATE
R2	sec	(B) µl	0 µl	0 µl	0 µl	
R3	0 sec	0 µl	0 µl	0 µl	NONE µl	
Photometry	Time	Cartridge configuration				
P1	-30 sec	Component	(A)	(A)	(A)	(A)
P2	300 sec	Number of tests	15	15	15	15
P3	sec	Well life hours	720	720	720	720
P4	sec	On board life 720 hours	Calibration 2160 hours			
Calculation						
Mode	END POINT	Measuring filter	540 nm	Blanking filter	700 nm	
P1	Time -30 sec	Dilution	-1.0	IOD	0	(Inactive)
P2	Time 300 sec	Dilution	1.0	FOD	0	(Inactive)
	1	2	3	4	5	6
	*	*	*	*	*	*
Decimal	2					
Units	g/L					
Interval	0.66 – 0.87					
Range	128					

Note: * entered by operator.