

# DADE DIMENSION PARAMETER SETTINGS

# G-cell

## UREA

Channel:	*	Name:	XURE	Mode:	ABSORBENCE	Std. curve:	LINEAR
<b>Delivery</b>							
	<b>Time</b>		<b>Component 1</b>	<b>Component 2</b>	<b>Component 3</b>	<b>Chase</b>	<b>Mix</b>
R1	-60.0 sec	(A)150 µl	0 µl	0 µl	0 µl	20 µl	MODERATE
S1	0.0 sec	2.5 µl	0 µl	0 µl	0 µl	20 µl	MODERATE
R2	180 sec	(B)50 µl	0 µl	0 µl	0 µl	20 µl	MODERATE
R3	0 sec	0 µl	0 µl	0 µl	0 µl	NONE µl	
<b>Photometry</b>							
	<b>Time</b>		<b>Cartridge configuration</b>				
P1	240 sec	Component	(A)	(A)	(A)	(A)	(B) (B)
P2	360 sec	Number of tests	15	15	15	15	30 30
P3	sec	Well life hours	720	720	720	720	720 720
P4	sec	On board life	720 hours	Calibration	2160 hours		
<b>Calculation</b>							
<b>Mode</b>	RATE	<b>Measuring filter</b>	340 nm	<b>Blanking filter</b>	405 nm		
P1	Time 240 sec	Dilution	1.0	IOD	0	(Inactive)	
P2	Time 360 sec	Dilution	1.0	FOD	0	(Inactive)	
	1	2	3	4	5	6	
	*	*	*	*	*	*	
<b>Decimal</b>	2						
<b>Units</b>	Mg/dL						
<b>Interval</b>	10-50						
<b>Assay Range</b>	258						