

Specialised Parameter Medicon Reagent Application sheet for Transasia Erba Analyzer Serises.

Application Sheet Code : RCD/ERBA/MED/001

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 |
|---------------------|------------------|-------------|------------|--------------------|-------------------|---------------------|--------------------|----------------|-------------------|-----------------------------------|
| Code | 3511-0234 | 1418-0510 | 1418-1042 | 1417-0527 | 1418-0274/0273 | 1417-0910 | 1418-0800 | 1417/1418-0120 | 1418-0740 | 1418-0280 |
| Method ID | ACP | ALD | AMO | CrP | FER | G6PD | IgE | IRON | TRF | UIBC |
| Product Name | Acid Phosphatase | Aldolase | Ammonia | C-reactive Protein | Ferritine | Glucose-6-Phosphate | Immunoglo-buline E | IRON | Transferrien | Unsaturated Iron Binding Capacity |
| Method | Fast-Red | UV Kinetics | Enzymatic | Latex Enhanced IT | Latex Enhanced IT | UV Kinetics | Latex Enhanced IT | Ferrozoin | Latex Enhanced IT | Nitrozo-PSAP |
| Test Detail | | | | | | | | | | |
| Test | ACP | ALD | AMO | CrP | FERR | G6PD | IgE | IRON | TRF | UIBC |
| Host Name | RCD | RCD | RCD | RCD | RCD | RCD | RCD | RCD | RCD | RCD |
| Report Name | ACP | ALD | AMO | CrP | FERR | G6PD | IgE | IRON | TRF | UIBC |
| Unit | U/L | U/L | µmol/L | mg/dL | ng/mL | U/L | IU/mL | µg/dL | mg/dL | µg/dL |
| Decimal Places | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 |
| Wavelength (nm) | 410 | 340 | 340 | 340 | 800 | 340 | 600 | 570 | 520 | 800 |
| Primary | ** | ** | 380 | 800 | ** | 380 | ** | 800 | 800 | ** |
| Secondary | | | | | | | | | | |
| Assay Type | Rate - A | Rate A | 2 - Point | 2 - Point | 2 - Point | Rate -A | 2 - Point | 2 - Point | 2 - Point | 2 - Point |
| Curve Type | Linear | Linear | Linear | Polynomial | Cubic Spline | Linear | Polynomial | Linear | Polynomial | Linear |
| *EM200 | | | | | | | | | | |
| M1 Start | 0 | 0 | 15 | 15 | 0 | 0 | 0 | 15 | 15 | 15 |
| M1 End | 0 | 0 | 16 | 16 | 0 | 0 | 0 | 16 | 16 | 16 |
| M2 Start | 9 | 13 | 35 | 35 | 19 | 21 | 19 | 29 | 35 | 35 |
| M2 End | 14 | 27 | 36 | 36 | 27 | 31 | 29 | 30 | 36 | 36 |
| *EM360 | | | | | | | | | | |
| M1 Start | 0 | 0 | 11 | 11 | 0 | 0 | 0 | 11 | 11 | 11 |
| M1 End | 0 | 0 | 12 | 12 | 0 | 0 | 0 | 12 | 12 | 12 |
| M2 Start | 16 | 24 | 50 | 50 | 17 | 20 | 16 | 37 | 50 | 50 |
| M2 End | 25 | 48 | 51 | 51 | 43 | 38 | 34 | 38 | 51 | 51 |
| *XL640 | | | | | | | | | | |
| M1 Start | 0 | 0 | 23 | 23 | 0 | 0 | 0 | 23 | 23 | 23 |
| M1 End | 0 | 0 | 24 | 24 | 0 | 0 | 0 | 24 | 24 | 24 |
| M2 Start | 18 | 27 | 63 | 63 | 30 | 33 | 28 | 51 | 63 | 63 |
| M2 End | 28 | 54 | 64 | 64 | 60 | 53 | 48 | 52 | 64 | 64 |
| Sample Replicates | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 |
| Standard Replicates | 3 | 3 | 3 | 3 | 3 | 3 | 3 | 3 | 3 | 3 |
| Control Replicates | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 |
| Control intervals | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 |
| Reaction Direction | Increasing | Decreasing | Decreasing | Increasing | Increasing | Increasing | Increasing | Increasing | Increasing | Decreasing |
| React Abs Limit | 2.5 | 0 | 0 | 2.5 | 2.5 | 2.5 | 2.5 | 2.5 | 2.5 | 0 |
| Prozone limit % | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
| Prozone Check | Lower | Lower | Lower | Lower | Lower | Lower | Lower | Lower | Lower | Lower |
| Delta Abs/Min | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
| Technical Minimum | 0.3 | 0.37 | 5.3 | 0.13 | 6 | 2.2 | 5.7 | 1.5 | 3 | 21 |
| Technical Maximum | 60 | 75 | 700 | 30 | 450 | 3000 | 1500 | 1000 | 600 | 400 |
| Y=aX+b | a= | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 | 1 |
| | b= | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
| Reagent Abs Min | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
| Reagent Abs Max | 1 | 0 | 0 | 1 | 1 | 1 | 1 | 1 | 1 | 0 |
| Test Volume | | | | | | | | | | |
| Sample Volumes | | | | | | | | | | |
| Normal | 40 | 16 | 60 | 23 | 24 | 12 | 6 | 40 | 3.6 | 9 |
| Increase | 40 | 32 | 60 | 23 | 48 | 12 | 12 | 40 | 3.6 | 18 |
| Decrease | 20 | 16 | 60 | 23 | 24 | 12 | 6 | 40 | 3.6 | 4.5 |
| Standard Volume | 40 | 16 | 60 | 23 | 24 | 12 | 6 | 40 | 3.6 | 9 |
| RGT-1 Volume | 200 | 200 | 180 | 150 | 160 | 300 | 150 | 180 | 180 | 186 |
| RGT-2 Volume | ** | ** | 60 | 60 | 160 | 25 | 48 | 180 | 180 | 30 |